

EPROM ERASURE KIT.


Kit includes Ballast Unit Starter, 12" 8W. 2537 Angst U.V. Tube, pair of bi-pin leads, Neon Indicator, safety Microswitch, on/off Switch and Circuit (Less Case).

Typical Erasure time is 3 to 5 minutes with the I.C. 1" to 2" from the Tube. As the usable area of the Tube is approx. 10" up to 8 or 10 I.Cs. may be treated simultaneously.

NOTE: The light output of this Tube is HIGHLY DANGEROUS to the eyes and every precaution must be taken to avoid exposure.

A suitable Box must be constructed and a Microswitch fitted so that the moment the lid is raised the lamp will be extinguished.

CIRCUIT DIAGRAM



Insert tube into bi-pin holders and twist to secure.

Although a Box with lifting lid is suitable, a more practical Box would be where the I.Cs. are first placed on a tray with open slot, and the tray slid into an opening in the Box with safety microswitch fitted at the far end so that the lamp does not switch on until the tray is fully inserted.

We must insist that under no circumstances may the Tube be switched on where it's light may be seen by humans or animals.

All Mail Orders—Callers—Ample Parking
Dept. PL, 57 BRIDGMAN ROAD,
CHISWICK, LONDON W4 5BB
Phone: 01-995 1560
Showroom open Mon—Fri

**SERVICE
TRADING CO.**

Personal callers only. Open Sat.
9 LITTLE NEWPORT STREET,
LONDON WC2H 7JJ
Phone: 01-437 0576